



Intravenous news: FDA report “A number of errors have been reported by Abacus software users as a result of ordering salt based parenteral nutrition ingredients on an ion based ordering template. Abacus TPN Calculation Software is designed and intended to allow the ordering of electrolytes in only one of two ways: as a salt (such as calcium gluconate 10%) or as an elemental ion (such as calcium). However, if a dosage is entered into the system based on one method, when the template is configured for the other method, a dosing error can occur.”

[Click here for the full story.](#)



Baxa Corporation Abacus Total Parenteral Nutrition (TPN) Calculation
Software: Class I Recall - Potential dosing errors | 2



Baxa Corporation Abacus Total Parenteral Nutrition (TPN) Calculation
Software: Class I Recall - Potential dosing errors | 3

