



FDA report “The Plum A+ and A+3 infusion pumps have an alarm that should sound when a therapy is interrupted. Some of the alarms may fail to sound in situations that should trigger it. It is possible for a long delay before a health care professional becomes aware of the need to restore therapy. For patients receiving critical intravenous medication, there is a risk of injury or death resulting from this prolonged interruption in therapy.”

Plum A+ and Plum A+3 infusion systems alarm failure <http://ctt.ec/l2Ud0+> @ivteam #ivteam

Click To Tweet

Full Alert

**Thank you to our partners for supporting IVTEAM**

