



Bayer has determined all Source Administration Sets used with the Medrad Intego PET Infusion System may produce a particulate matter in the medicine vials” FDA (2017).

FDA (2017) report “Bayer has determined all Source Administration Sets used with the Medrad Intego PET Infusion System may produce a particulate matter in the medicine vials. The particulates may be created when the tip of the needle pushes through the rubber top of the vial. If this occurs, the particulate matter could enter into the patient and cause serious adverse health consequences including infection, damage of tissue, and death.”

ReTweet if useful... Recall of Medrad Intego PET Infusion System Source Administration Sets by Bayer <https://ctt.ec/9RrE0+> @ivteam #ivteam

Click To Tweet

Full Alert

Thank you to our partners for supporting IVTEAM

