

ChloroPrep® venous access and maintenance range

2% Chlorhexidine Gluconate w/v and 70% Isopropyl Alcohol v/v
Licensed, single use, sterile, ANTT™, latex free, skin antiseptics system

Pinch



Holding the barrel of the applicator, pinch the levers together, being careful not to touch the sponge.

Pinch the lever to release the solution. You will hear a 'pop' as the ampoule breaks.

Apply



Starting at the incision site, gently press the applicator against the skin until the solution soaks the sponge.

Apply using repeated up and down, back and forth strokes for at least 30 seconds, before working outwards to the periphery.

Dry



Leave the area to air **dry completely** before applying sterile drapes. Do not blot or wipe away. Discard the applicator after a single use.

Important safety point

Do not drape or use ignition source until the solution has completely dried.

3ml applicator

Coverage area
15 cm x 15 cm



For procedures such as:

- Central venous catheter insertions
- Insertion site maintenance
- Midline catheter insertions
- Minor procedures

10.5ml applicator

Coverage area
25 cm x 30 cm



For procedures such as:

- Surgical site antisepsis
- Tunelled central venous catheter insertions
- Interventional radiology
- Implantable device placements