

Chloraprep® SEPP .67ml cutaneous solution

2% Chlorhexidine Gluconate w/v and 70% Isopropyl Alcohol v/v

Licensed, single use, sterile, ANTT™, latex free, skin antiseptics system

Pinch



Hold the applicator between thumb and forefinger, being careful not to touch the tip. **Pinch** the sides of the applicator once (**in the middle of the ampoule**) until you feel the inner ampoule break.

Do not keep 'pinching'
Do not bend or snap the applicator

Apply



Let the solution flow naturally through the non-linting tip.

Gently press the applicator against the skin and apply the antiseptic using repeated up and down, back and forth strokes for about 30 seconds.

Dry



Leave for approximately 30 seconds, allowing the area to air **dry completely** before proceeding.

Discard the applicator after a single use.

SEPP applicator

Coverage area
5 cm x 8 cm



For procedures such as:
Peripheral cannulation
Surgical pin site cleansing
Subcutaneous device insertions

For customer services and all other enquiries: Please telephone: 0800 043 7546

Or visit: www.chloraprep.co.uk

CHL065a Sept 2008

Chloraprep®
Patient Preoperative Skin Preparation
2% Chlorhexidine Gluconate (CHG) w/v
& 70% Isopropyl Alcohol (IPA) v/v