

Think Positive

BD™ Saline

is changing name to

BD PosiFlush™

BD PosiFlush™ is a range of ready to use sterile pre-filled syringes containing sodium chloride.

They are specifically designed for flushing in-situ vascular access devices in order to maintain patency and reduce the risk of contamination.



All the **Positive** benefits for the safe and optimal **Flush** of vascular access devices remain:

- Clear labelling for easy identification and verification
- Minimise the use of potentially contaminated flushing solutions
- Designed to prevent solution from entering a non-sterile area of the syringe
- Available for both sterile field (XS) and non-sterile field (SP) applications
- Special design to eliminate syringe-induced blood reflux
- Special syringe design reduces storage and disposal cost and minimises environmental waste
- Special tip cap (with screw thread) minimises contamination risk
- Reduces the number of steps and time to prepare a flush syringe, creating time savings and efficiency

All catalogue numbers remain the same:

Product Description	Brand	BD Code	Units
<i>Syringe pre-filled flushing device sterile content and fluid pathway</i>			
Saline 3 mL fill volume	BD PosiFlush™ SP Syringe	306573	Box of 30
Saline 5 mL fill volume	BD PosiFlush™ SP Syringe	306574	Box of 30
Saline 10 mL fill volume	BD PosiFlush™ SP Syringe	306575	Box of 30
<i>Syringe pre-filled flushing device externally sterile</i>			
Saline 3 mL fill volume	BD PosiFlush™ XS Syringe	306570	Box of 30
Saline 5 mL fill volume	BD PosiFlush™ XS Syringe	306571	Box of 30
Saline 10 mL fill volume	BD PosiFlush™ XS Syringe	306572	Box of 30

